

MDE

MEDICAL DEVELOPMENT & ENGINEERING

HANDYCAST SAW



TECHNICAL SPECIFICATIONS

Supply	90 – 270 VAC, 50 – 60 Hz
Power	100 Watts
Weight	750 gm
Length	220mm
Working Diameter	47mm
Cable Length	Mains to power unit: 1 m
Cable Length	Power unit to Saw : 2m

ORDERING INFORMATION

➤ Handycast 2010 Saw	LH 1000 0000
➤ Blades – SC65	LH 1000 2001
➤ Blades – SC65C	LH 1000 2002
➤ Blades – SC50	LH 1000 2003
➤ Blades – SC50C	LH 1000 2004
➤ Allen Key	LH 1000 3000

SPECIAL FEATURES AND BENEFITS

- Ergonomically contoured saw designed to fit comfortably in hand while providing ideal cutting position
- Specially designed mechanism results in low noise level and reduced vibration
- Low weight (750 gm) & proper balancing ensures increased manoeverability & reduces wrist fatigue
- Maintenance free system ensures excellent reliability
- Completely sealed system for extensive use
- Incorporates state of the art brushless motors for unmatched performance & speed
- High performance blades for all types of application



SC65



SC65C



SC50



SC50C

MDE

MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026
 Tel : +91-20-27121512, 30720117 Fax. : +91-20-27121512
 Email : mdeindia@gmail.com Web: www.mdeindia.net

