

MICRO- DRILL MULTITASK

MDE

MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026

Tel : +91-20-27121512, 30720117, 30690147/8

Email : mdeindia@gmail.com Website : www.mdeindia.net



Industry's best and only completely Modular micro Power Orthopedic System.

Quick lock system enables attachment switching during a procedure without disconnecting hoses/cables.



- Powered by Micro Electric Console.
- Incorporates Powerful Brushless DC Motor for max torque in its class
- Cables & Connectors are Autoclavable
- Quick Lock Mechanism
- Can drive Electric Dermatome



MDE

MDE

MDE

MDE

MDE

MDE

MICRO-DRILL
MULTITASK

Micro Multitask System



Control Unit
Part No : LM 10000000



Hand Piece
Part No : LM 10020000



Reciprocating Saw
Part No : LM 10060000

Common Attachments



High Speed Drill attachment
Part No : LM 10030005



Drill Adapter
Part No : LM 10030006



Drilling attachment with
Mini Jacob's Chuck
Part No : LM 10030001



Foot Switch
Part No : LM 10010000




Saw Blade:
15x4x0.4mm
LM 10060001



Saw Blade:
25x4x0.4mm
LM 10060002



Micro Sagittal Saw
Part No : LM 10040000



Saw Blade:
15x6x0.4mm
LM 10040001



Saw Blade:
15x9x0.4mm
LM 10040002



Sagittal Saw(HALL Type)
Part No :LM 10050000



Saw Blade:
15x4.5x0.4mm
LM 10050001



Saw Blade:
25x9x0.4mm
LM 10050002



Saw Blade:
35x6x0.4mm
LM 10050003